

# Writing instruments made of antibacterial plastic.




ISO 22196 certified,  
hygienic and offering product protection for life.

**Be safe.**

**prodir®**

Antibacterial, hygienic surfaces are nowadays particularly important. That is why we have developed a whole portfolio for you right away. With new models, exciting new combinations and options that allow you to communicate what really counts:



-  +  **antibacterial**
- + DS3, DS5, DS8
- + 100% recycled
- + from 500 pieces
- + white

The new DS Regeneration Pens: sustainably ecological and sustainably antibacterial.

-  **antibacterial**
- + DS3
- + from 500 pieces
- + 30% recycled
- + white

The DS3 special edition: the icon, simple and convincing, in polished or matt.


-  **antibacterial**
- + all DS and QS models
- + from 5'000 or 10'000 pieces
- + 30% recycled
- + standard or Pantone colours

Your writing instrument, just as you like it, with an antibacterial surface.


The plastic surfaces of our writing instruments are protected by Biomaster antimicrobial/antibacterial technology. Biomaster is the world's most trusted additive, tested and certified to ISO 22196 providing durable and effective surface protection 24/7 for the lifetime of the product.



## How does Biomaster work?



**1** Biomaster binds the cell walls disrupting growth.



### How effective is Biomaster?

Very effective. Biomaster is proven to inhibit the growth of bacteria by up to 99.99%.

### Is Biomaster effective against antibiotic-resistant bacteria?

Yes. Biomaster has been proven to be effective to antibiotic-resistant bacteria.

### Is Biomaster safe?

Yes. It is based on silver ion technology, recognised for centuries with no harmful effects. Biomaster is used in medical, food and water applications.

### Does Biomaster use nano-silver technology?

No. No nano technology is used due to ongoing safety concerns.

### Has Biomaster been tested?

Yes, repeatedly. Biomaster is tested to ISO standards and is subject to continuous quality control tests and environmental trials.

### Why is Biomaster better than other silver-based antimicrobial additives?

Biomaster pioneered the use of antimicrobial additives and remains the recognised leader and most trusted supplier of antimicrobial technology for polymers, textiles, paper, paints and coatings.

### How long is Biomaster effective?

Biomaster is effective for the whole lifetime of the pen.


Please note, that Biomaster is tested to ISO 22196 providing efficacy against a wide range of microbes including bacteria. Efficacy against ISO 21702 for viruses is ongoing.

### Does Biomaster affect a product in any way?


No. You can't see, smell or even taste Biomaster.

More about Biomaster [here!](#)

**2** The Biomaster ions interfere with the enzymes production stopping the cell producing energy.



**3** Biomaster interrupts the cell's DNA preventing replication.



Effective and save.



# New Antibacterial DS3 DS5 DS8 Regeneration Pens



The unique combination of sustainability and hygiene with antibacterial body made of recycled plastic. Available in the standard version in fine white matt from 500 pieces.



Colours

Button (DS8)

Regenerated antibacterial



J02

Clip (DS8)

Regenerated antibacterial



J02


Cap and barrel

Regenerated antibacterial




N02

DS3 Print areas



DS5 Print areas



DS8 Print areas



Minimum order quantity  
500 pens

 For PMS with antibacterial protection starting from 5.000 pieces

	mm	Inch	Print
<b>C</b> Back cap	25 × 4	0,98 × 0,16	1 colour
<b>B</b> Clip	40 × 8	1,57 × 0,31	6 colours
<b>D</b> Body	50 × 25	1,97 × 0,98	4 colours

	mm	Inch	Print
<b>C</b> Back cap	23 × 5	0,90 × 0,20	1 colour
<b>B</b> Clip	40 × 7	1,57 × 0,28	6 colours
<b>D</b> Body	35 × 30	1,38 × 1,18	4 colours

	mm	Inch	Print
<b>B</b> Clip	36 × 7	1,41 × 0,28	6 colours
<b>C</b> Body	36 × 5	1,41 × 0,20	2 colours
<b>C1</b> Side body	36 × 5	1,41 × 0,20	2 colours
<b>C2</b> Side body	36 × 5	1,41 × 0,20	2 colours



DS5 TNN

Cap N02 Antibacterial  
 Body N02 Antibacterial  
 Nose cone N02 Antibacterial

DS3 TNN

Cap N02 Antibacterial  
 Body N02 Antibacterial  
 Nose cone N02 Antibacterial

DS5 TNN

Cap N02 Antibacterial  
 Body N02 Antibacterial  
 Nose cone N02 Antibacterial

DS8 PNN-J

Cap N02 Antibacterial  
 Barrel N02 Antibacterial  
 Clip J02 Antibacterial  
 Button J02 Antibacterial



DS8 PNN-J


Cap Matt PMS Antibacterial  
 Barrel Matt PMS Antibacterial  
 Clip Polished PMS Antibacterial  
 Polished PMS Antibacterial



# New Antibacterial DS3



The icon, simple, clear and timeless. Available in the antibacterial standard versions in matt or polished white already from 500 pieces.



## Colours

### Cap, body and nose cone

Matt  
Antibacterial




M02

Polished  
Antibacterial



P02

## DS3 Print areas



Minimum order quantity  
500 pens



For PMS with antibacterial  
protection starting from 5.000 pieces

	mm	Inch	Print
<b>B</b> Clip	40 × 8	1,57 × 0,31	6 colours
<b>C</b> Back cap	25 × 4	0,98 × 0,16	1 colour
<b>D</b> Body	50 × 25	1,97 × 0,98	4 colours



DS3 TMM

Cap M02 Antibacterial  
Body M02 Antibacterial  
Nose cone M02 Antibacterial

DS3 TMM

Cap M02 Antibacterial  
Body M02 Antibacterial  
Nose cone M02 Antibacterial

DS3 TPP

Cap P02 Antibacterial  
Body P02 Antibacterial  
Nose cone P02 Antibacterial

DS3 TPP

Cap P02 Antibacterial  
Body P02 Antibacterial  
Nose cone P02 Antibacterial



DS3 TPP

Cap Polished PMS Antibacterial  
Body Polished PMS Antibacterial  
Nose cone Polished PMS Antibacterial



New Antibacterial

# All DS models All QS models



All models are available with antibacterial finish in both standard and special colours from a minimum order quantity as the one of the special colours of each model.\*



\*Please refer to the next page for more information.



# Good to know

.....

Writing instruments with antibacterial surface protection can be recognised by the imprint or the relief „antibacterial“ on the surface.

Execeptions are the models QS40 Air and QS01 which have a single print area on the clip. In this case, to identify the pen as antibacterial, we recommend the use of packaging or a simple cardboard case.

.....

Antibacterial surface protection is not available for DS3 Biotic and QS Stone models, for varnished or Soft touch surfaces, and for metal components such as clips, nose cones and buttons.

.....

All writing instruments are equipped ex works with high-quality, low-emission Floating Ball® 1.0 mm Lead Free refills. Other refills from the 2020 catalogue are optionally available.

.....

The casing of all our writing instruments is made from at least 30% recycled plastic.

For our Regeneration Pens the percentage is 100%.



# Where we are

## Headquarters

Pagani Pens SA  
BU Pens - Prodir  
Via Serta 22  
CH-6814 Lamone  
Switzerland

Tel +41-91-935 55 55  
Fax +41-91-935 55 35  
sales@paganipens.com

## Prodir Germany

Prodir GmbH  
An der Molkerei 8  
D-56288 Kastellaun  
Germany

Tel +49-6762-4069-0  
Fax +49-6762-4069-69  
sales.de@paganipens.com

## Prodir Austria

Prodir Austria GmbH  
Hegelgasse 8 / 10  
AT-1010 Wien  
Austria

Tel +43-7205 69480  
Fax +43-2622-40301-15  
austria@paganipens.at

## Prodir Benelux

Prodir BV  
Fokkerstraat 28  
NL-3905 KV Veenendaal  
Netherlands

Tel +31-318-74 42 70  
Fax +31-318-51 48 44  
sales.nl@paganipens.com

## Prodir Scandinavia

Prodir Aps Scandinavia  
Ambolten 8, 1.th.  
DK-2970 Hørsholm  
Denmark

Tel +45-88-44 12 50  
Fax +45-48-24 57 90  
sales@paganipens.dk

## Prodir Great Britain

Prodir Ltd  
5 South Fens Enterprise Park  
Fenton Way  
Chatteris, Cambridgeshire  
PE16 6WA  
UK

Tel +44-1354-708590  
Fax +44-1354-776490  
sales@paganipens.co.uk

## Prodir France

Prodir France SAS  
5 rue Charles de Gaulle  
94140 Alfortville  
France

Tel +33-1-79 35 74 30  
Fax +33-1-42 33 00 40  
sales@paganipens.fr

## Prodir Spain

Prodir Iberia SA  
C/Copenhagen, 12  
Edificio Tifán Oficina 213  
Polígono Industrial Európolis  
ES - 28232 Las Rozas  
Spain

Tel +34-91-142 08 40  
Fax +34-91-522 14 26  
sales@paganipens.es

## Prodir Russia

Prodir SA  
Representative Office  
Smolenskaya square 3  
Moscow – 121099  
Russia

Tel +7-495-782-1915  
nushin.loyko@paganipens.ru

